



[-> to the Product](#)

PEPLYN PLUS filter elements are made entirely of polypropylene and are used primarily in the beverage industry or also in applications for pre-stabilization of process liquids. The filter elements with absolute separation rate represent an optimal and cost-effective solution for treating a wide variety of beverage products prior to intermediate storage. The continuously graded fiber density represents a major advantage in filtration, as it enables increasingly finer particles to be retained in the filter depth. In combination with an optimum pleat density, the PEPLYN-PLUS offers excellent separation performance and a very long service life.

Main features


- Absolute separation efficiency from 0.6 µm to 100 µm
- Pleated filter element with graduated fiber density
- Materials conform to the relevant requirements of FDA 21 CFR Part 177, current 1935/2004 and current USP Plastics Class VI - 121 °C
- High flow rate and long service life
- Can also be used under aggressive process conditions

Traceability

- Individual serial number, LOT number and product code for each element

Sterilization

- Steam sterilizable and autoclavable up to 135 °C
- Hot water sterilization with up to 90 °C
- Compatible with a wide range of chemicals



Technical data


Material:	Polypropylene
Insert end caps:	316L stainless steel *not available for adapter B
Micron rating:	0,6 µm to 100 µm
Gasket:	Silicone (Standard), others on demand
Element length:	10" - 40"
Filter surface:	0,65 m ² (per 10"- Filter element)
Max. working temperature:	70°C

Separation efficiency depending on micron rating

efficiency	> 99,99%	99,98%	99,90%	99%	90%
β	10000	5000	1000	100	10
0,6	0,6	0,57	0,54	0,32	0,2
1,0	1,0	0,95	0,9	0,7	0,5
1,5	1,5	1,4	1,1	0,8	0,6
3	3	2,8	1,8	1,0	0,7
5	5	4,7	4,5	3,5	1,0
7	7	6,7	6,3	4,5	2,5
10	10	8	7	4,8	2,8
15	15	12	10	7,2	4,5
20	20	16	14	10	6
25	25	20	17	12	7

Short time temperature for CIP

Temperature [°C]	Differential pressure forward [bar]
20	5,0
40	4,0
60	3,0
80	2,0
90	1,0
>100 (steam)	0,3



Order information

ZCPP	[]	-	[]	-	[]	-	[]	-	[]
Length	Micron rating	Adapter	Variant	Gasket					
1 10"	.60 0,6 µm 015 15 µm	B* dh DOE	S* steam-sterilizable	S* Silicone					
2 20"	1.0 1,0 µm 020 20 µm	C BF/226/Bajonett							
3 30"	1.5 1,5 µm 025 25 µm	D Fin/222							
4 40"	003 3,0 µm 040 40 µm	N Style 216/Code6							
	005 5 mm 055 55 mm								
	007 7 µm 075 75 µm								
	010 10 µm 100 100 µm								

Order example: ZCPP3-1.0CS